



AEGIS LIFESCIENCES

PROTECTION THROUGH PERFECTION

CE
2460

FDA
REGISTERED FACILITY

ISO
13485

MDSAP

WE
MATTER
WHEN
HAEMOSTATS
MATTER




Aegis Lifesciences, established in the year 2008, is world's leading manufacturer & exporter of all major types of Haemostats. It has gained **international certificates of CE 2460, MDSAP & ISO 13485 Quality Management System. Exporting to 100+ countries** including highly regulated regions like Europe (CE), USA, Canada, Brazil (ANVISA), UK, Russia, Middle East etc.

MISSION

To improve the lives of patients around the world by developing and manufacturing innovative, high-quality & cost-effective medical devices.

VISION

To be the leading global provider of hemostatic products & technologies.





SURGISPON®
ABSORBABLE HAEMOSTATIC GELATIN SPONGE



Surgi-ORC®
OXIDIZED REGENERATED CELLULOSE HAEMOSTAT



FLOGEL®
FLOWABLE GELATIN HAEMOSTATIC MATRIX KIT



SURGI-PVA™
POLYVINYL ALCOHOL MEDICAL SPONGE DRESSINGS



DENTOSPON®
ABSORBABLE GELATIN SPONGE WITH
5% COLLOIDAL SILVER



HAEMOVET®
ABSORBABLE HAEMOSTATS FOR VETERINARY USE

SURGISPON®

ABSORBABLE HAEMOSTATIC GELATIN
SPONGES & POWDER



Scan For Detailed
Information

CHARACTERISTICS

- Absorbs up to 45 times its own weight
- pH-neutral, ideal carrier for drugs
- Completely absorbed within 4 weeks
- Can be used dry or saturated with sterile saline solution
- Non-toxic, non-allergenic, non-immunogenic, and non-pyrogenic
- Uniform porosity guarantees a favourable haemostasis
- No residue & no encapsulation post complete absorption
- Easy to cut according to need
- Extremely cost-effective
- No special storage conditions required

**AVAILABLE IN 30+
CUSTOMIZED SIZES
& POWDER FORM**



SURGISPON[®]

ANAL TAMPON



Scan For Detailed
Information

CHARACTERISTICS

- Reduces post-operative pain
- No need for manual removal
- Can reduce up to 4 days of hospital stay
- Easy expelled at first bowel movement
- Liquefies rapidly and is spontaneously expelled within 1 to 2 days
- Painless and easy tampon discharge
- Provided with an opening into which flatus drain can be introduced
- Best suitable for post-operative haemostasis
- Well-adhering ability to the bleeding site

AVAILABLE SIZES

80xØ30 mm

SURGISPON®

BELLOW POWDER (DRY USAGE)



Scan For Detailed
Information

CHARACTERISTICS

- Stops bleeding effectively with significantly faster time
- For large and irregular surfaces
- Significantly lower chances of rebleeding
- Does not leave voids between the clot material and tissue surfaces
- Provides better mechanical integrity to the clot by enhancing the fibrin matrix within clots
- Better blood absorption capacity
- No risk of encapsulation and overdosage
- For use in open and laparoscopic surgeries
- Can be easily removed or irrigated once haemostasis is achieved

AVAILABLE SIZES

1 g, 2 g, 3 g, 5 g
14 cm & 38 cm laparoscopic
applicators



SURGISPON[®]

CONTAINER POWDER (PUTTY)




Scan For Detailed
Information

CHARACTERISTICS

- Can be easily formed into any shape, according to surgical site
- Applied dry in powder form or as a paste mixed with saline solution
- Used in sternotomy and craniotomy to control bleeding in bones, cavities, and fissures
- Can be applied by hand or with syringe to areas that are hard to reach
- Does not cause bone dehiscence like conventional bone waxes
- Fast & proven haemostasis for those hard-to-reach areas & oozing spongy bone

AVAILABLE SIZES

1 g, 3 g & 5 g



SURGISPON®

DENTAL RANGE



Scan For Detailed
Information

CHARACTERISTICS

- Can be used as a carrier matrix along with drugs in bone regeneration
- User-friendly: easy to trim and manipulate
- Allows antibiotic addition for high-risk patients
- Mitigation of alveolar ridge resorption
- Saves valuable surgical and anaesthetic time
- Used as a scaffold for slow release of growth factors in healing of extraction sockets
- Supports stable blood clot formation

AVAILABLE SIZES

10x10x10 mm, 12xØ8 mm,
20xØ10 mm & 20x20x07 mm
1 g Bellow powder & 1 g Putty powder

DENTOSPON[®]

ABSORBABLE GELATIN SPONGE WITH
5% COLLOIDAL SILVER



Scan For Detailed
Information

CHARACTERISTICS

- Smooth & complete resorption within 4 weeks
- Long lasting antimicrobial depot effect
- Supports stabilization of the blood clot
- Complication-free healing
- It is both haemostatic and bactericidal
- Finely dispersed colloidal silver provides active surface for continuous release of silver ions
- Colloid silver Ph. Eur. 0.5 mg
- Reliable post-extraction treatment
- Effective against all micro-organisms found in oral cavity, particularly against gram-negative organisms

AVAILABLE SIZES

10x10x10 mm, 14x7x7 mm

Surgi-ORC®

OXIDIZED REGENERATED CELLULOSE HAEMOSTAT
ORIGINAL/STANDARD



CHARACTERISTICS

- Fine weave of Oxidized Regenerated Cellulose
- Sheer weave allows visibility of the bleeding site
- Complete bio-absorption within 7-14 days
- pH<4, provides bactericidal properties
- Can be cut without fraying
- High flexibility and drapeability
- Suitable for both open and MIS surgeries
- Light knit formation covers uneven tissue
- Safe haemostasis with effective bactericidal properties

AVAILABLE SIZES

2x14 in, 4x8 in, 2x3 in,
0.5x2 in



Scan For Detailed
Information

Surgi-ORC®

OXIDIZED REGENERATED CELLULOSE HAEMOSTAT FIBRIL



CHARACTERISTICS

- Tufts of soft, lightweight ORC
- 7 layers of fibrillose fabric
- Provides high flexibility and conformability to surgical sites
- Soft, layered material conforms and “melts in” to bleeding tissue
- Lightweight layers can be pulled apart for precise placement around bleeding site
- Convenient for hard-to-reach site or irregularly shaped bleeding site
- Layers can be “peeled” off in desired amounts

AVAILABLE SIZES

1x2 in, 2x4 in, 4x4 in



Scan For Detailed
Information

Surgi-ORC®

OXIDIZED REGENERATED CELLULOSE HAEMOSTAT
NON-WOVEN/SNOW



CHARACTERISTICS

- Non-woven design increases contact surface
- Needle punched with interlocking fibres
- Easy use in open and MIS procedures without sticking to the instruments
- Faster Haemostasis than original variant
- Minimal to mild bleeding from specific or widespread surgical regions
- Retains its shape to offer conformability, handling & adherence to tissue, especially in laparoscopic surgery

AVAILABLE SIZES


1x2 in, 2x4 in, 4x4 in



Scan For Detailed
Information

Surgi-ORC®

OXIDIZED REGENERATED CELLULOSE HAEMOSTAT KNIT



Scan For Detailed
Information

CHARACTERISTICS

- Denser weave of ORC for heavier bleeding
- Soft, pliable weave designed to hold a suture
- Supports stabilization of the blood clot
- Ideal web structure where wrapping and suturing is required
- Dense ORC material provides high tensile strength

AVAILABLE SIZES

2x3 in, 3x4 in, 6x9 in,
1x1 in, 1x3 in, 1x3.5 in

Surgi-ORC®

OXIDIZED REGENERATED CELLULOSE HAEMOSTAT POWDER



Scan For Detailed
Information

CHARACTERISTICS

- Natural, plant-based absorbable haemostat
- Dual-particle size composition – “Finer ones” & “Coarser ones”
- Faster stopping of continuous oozing on broad surfaces
- Forms a more adherent and durable clot
- Maintains haemostasis even under irrigation
- Low pH confers a bactericidal effect, enhancing patient safety
- Achieves haemostasis in < 1 minute
- Free from animal-origin substances or collagen additives
- For use in open & laparoscopic surgeries

AVAILABLE SIZES

1 g, 2 g, 3 g, 4 g, 5 g
14 cm & 38 cm laparoscopic applicators



FLOGEL®

FLOWABLE GELATIN HAEMOSTATIC MATRIX KIT



CHARACTERISTICS

- To stop bleeding in tight & irregular spaces, in difficult-to-reach surgical site
- Envelops bleeding site to stop bleeding faster
- Maintains uniform viscosity from beginning to end
- Can be mixed with either saline or thrombin
- Can be used up to 8 hours after mixing with the sterile saline or thrombin
- Uniform and consistent flowability for precise performance
- Absorbed completely in 4 weeks
- Less stringent storage requirements compared to competition

AVAILABLE SIZES

5 ml, 6 ml, 8 ml & 10 ml



Scan For Detailed
Information



SURGISPON®

ENT RANGE



Scan For Detailed
Information

CHARACTERISTICS

- Supports vital tissue structures
- Prevents adhesions and Minimizes edema
- Firm, higher density than standard cubes
- Stops bleeding faster and effectively
- Avoids risk of tissue ingrowth & bone formation
- Anti-adherence to the packing material, gloves & instruments
- Separates tissues & maintains moisture
- Drains prior to first post-op appointment for clear progress evaluation
- Painless clearance through natural mucosal flow

AVAILABLE SIZES

80x50x10 mm, 10x10x10 mm,
10x10x10 mm HD, 20x20x07 mm,
80x15x07 mm, 80xØ07 mm,
1 g Bellow powder & 1 g Putty powder

SURGI-PVA™

NASAL PACKS



Scan For Detailed Information

CHARACTERISTICS

- Effective and ideal solution to assist in management of fluids
- Expands when moistened to assist in fluid management
- Air-foaming technology, prevents bacterial growth, lowers infection & reduces odor
- Soft, supportive, anti-adhesive and low-irritation
- 3D open cell interconnected structure gives strength & softness to the sponge for atraumatic removal
- Provides optimal pressure during epistaxis treatment
- Breathable airway, effectively reduce patient discomfort
- Available with Thread and/or Airway Tube

AVAILABLE SIZES

Standard Nasal Packing - 45x20x15 mm & 80x20x15 mm
 Anatomical Nasal Packing - 80x30x15 mm
 Epistaxis Nasal Packing - 100x25x15 mm
 Sinus Nasal Packing - 35x12x6 mm & 35x12x9 mm

SURGISPON®

NASAL/DERMAL WOUND CARE




Scan For Detailed
Information

CHARACTERISTICS

- Hypoallergenic, sterile, natural and safe
- Stops bleeding effectively
- Softness allows for painless removal
- Strong ability to absorb blood and exudate
- pH-neutral - ideal carrier for antibiotics & other medicinal drugs
- For persistent nasal and dermal bleedings resulting from injuries, accidents, allergies etc.
- Insoluble in water and fully absorbable

AVAILABLE SIZE

40x10x10 mm



SURGISPON®

COMPRESSED NASAL STRIP




Scan For Detailed
Information

CHARACTERISTICS

- Absorbs up to 45 times its own weight
- pH-neutral, inert & an ideal drug carrier (antibiotics, thrombin etc.)
- Conforms to nasal cavity for easy placement & better visibility
- Expands on fluid absorption to provide effective tamponade effect
- Used as a topical haemostat and dressing
- Serves as an intranasal splint/pack
- Provides support and scaffold for tissue regeneration
- Helps prevent adhesions between the septum & the nasal cavity

AVAILABLE SIZES

80x15x02 mm, 40x10x02 mm



SURGISPON®

HAEMOSTATIC PATCH



Scan For Detailed
Information

CHARACTERISTICS

- Designed to minimize bleeding in the sutured area
- Improved suture retention, low thrombogenicity, and resistance to aneurysmal dilation
- Durable: Non-fraying with multi-directional strength
- Ease of Handling: Can easily be trimmed to the desired size, and accommodated by cardiovascular anatomy
- Compatible with antibiotics due to neutral pH level
- Available in powder putty variant to control cancellous bone bleeding

AVAILABLE SIZES

40x60x05 mm, 70x100x05 mm,
60x120x05 mm, 80x160x05 mm
120x120x05 mm, 150x200x05 mm,
125x80x05 mm, 250x80x05 mm

SURGISPON® DIAL

FOR POST-HAEMODIALYSIS
BLEEDING CONTROL



Scan For Detailed
Information

CHARACTERISTICS

- Designed for post-haemodialysis patients to stop bleeding at the puncture site with light pressure
- Provides mechanical matrix promoting platelet aggregation & coagulation, facilitating clot formation
- Highly Porous Structure absorbs up to 45 times its own weight
- Fast & effective haemostasis, minimizing recovery time
- No reported side effects
- Leaves no residue or encapsulation behind
- Easy application and painless removal
- Can be used dry & in combination with drugs/antibiotics

AVAILABLE SIZES

30x30x10 mm, 20x20x07 mm

HOW TO USE

Step: 1



Step: 2



Step: 3



Step: 4



SURGI-PVA™

EYE SPEAR



Scan For Detailed
Information

CHARACTERISTICS

- Constructed using highly absorbent, soft & highly retentive unique Polyvinyl Alcohol (PVA) Sponge
- Biocompatible, clean, lint-free, and fibre-free
- Made through Air Foam Technology
- Not susceptible to bacterial growth
- Excellent tensile strength
- A widened, ergonomic grip provides better control
- Fast wicking & fluid control
- Soaks fluid > 8 times of its dry weight
- Lint-free material ensures a clean surgical field

AVAILABLE SIZE

17x8x2 mm

* Each Individual foil containing 5 sterile pieces

HAEMOVET®

ABSORBABLE HAEMOSTATS FOR
VETERINARY USE



Scan For Detailed
Information

CHARACTERISTICS

- Provides effective, quick & easy haemostasis
- It's porosity enables it to absorb up to 45 times
- Biologically absorbed within 4 weeks
- Saves valuable surgical & anaesthesia time
- Used dry or saturated with saline
- Can be cut to the desired size & shape
- Allows a dry, clearly visible surgical field
- Improves recovery time for the patient

AVAILABLE SIZES

10x10x10 mm, 60x20x07 mm,
80x50x10 mm, 125x80x10 mm
1/3/5 g Putty Powder,
1/2/3/5 g Bellow Powder,
With 5% colloidal silver-
10x10x10 mm, 14x07x07 mm,
20x20x20 mm





Customised size products to take care of almost all types of Surgeries and convenience to surgeons



Only facility in the world for exclusive manufacturing of absorbable haemostatic products, with largest production

30%+ Annual growth



Most cost-effective solutions to the ailing humanities

100+ Country presence across the globe

CORPORATE OFFICE

A-504, Swati Trinity, Nr. Shantipura Circle, SP Ring Road,
Ahmedabad, Gujarat-380058, INDIA.

MANUFACTURING SITE

Aegis Lifesciences Private Limited
215/216, Mahagujarat Industrial Estate, Ahmedabad, Gujarat-382213, India.

+91 9512903131 +91 9512883131

marketing1@aegis-lifesciences.com / aegis.lifesciences@gmail.com



WEBSITE



YOUTUBE

AL/CAT/GEN/01/R0 | Date:05/2026

Important: The selection and use of the product should be based on the surgeon's clinical expertise and the specific requirements of the surgical site.
Please refer to the Instructions for Use (IFU) for complete information on instructions, contraindications, warnings, and precautions.